



TCF-19 Polyclonal Antibody

Cat No: HR1AP4919

For research use only

Overview

Product Name	TCF-19 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GenID?Human?	TCF19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TCF19; SC1; Transcription factor 19; TCF-19; Transcription factor SC1
Background	transcription factor 19(TCF19) <i>Homo sapiens</i> This gene encodes a protein that contains a PHD-type zinc finger domain and likely functions as a transcription factor. The encoded protein plays a role proliferation and apoptosis of pancreatic beta cells. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],