



## Sox-17 Polyclonal Antibody

Cat No: HR1AP4896

For research use only

### Overview

Product Name	Sox-17 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GenID?Human?	SOX17
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SOX17; Transcription factor SOX-17
Background	SRY-box 17(SOX17) <i>Homo sapiens</i> This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. [provided by RefSeq, Jul 2008],