



## EF-1? Polyclonal Antibody

Cat No: HR1AP4870

For research use only

### Overview

Product Name	EF-1? Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GenID?Human?	EEF1B2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EEF1B2; EEF1B; EF1B; Elongation factor 1-beta; EF-1-beta
Background	eukaryotic translation elongation factor 1 beta 2(EEF1B2) Homo sapiens This gene encodes a translation elongation factor. The protein is a guanine nucleotide exchange factor involved in the transfer of aminoacylated tRNAs to the ribosome. Alternative splicing results in three transcript variants which differ only in the 5' UTR. [provided by RefSeq, Jul 2008],