

Histone H2B Mouse mAb

Cat No: HR1AM2274

For research use only

Overview

Product Name	Histone H2B Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-3,000 IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	14kDa
GeneID?Human?	8349
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	H2B GL105 antibody, H2B.1 antibody, Histone H2B type 1B antibody, HIST1H2BB antibody
Background	Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells.