



ZNF785 Polyclonal Antibody

Cat No: HR1AP5981

For research use only

Overview

Product Name	ZNF785 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	ZNF785
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF785; Zinc finger protein 785
Background	function:May be involved in transcriptional regulation.,sequence caution:The predicted gene has been split into 2 genes: ZNF688 and ZNF785.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 2 C2H2-type zinc fingers.,similarity:Contains 7 C2H2-type zinc fingers.,