



ZNF771 Polyclonal Antibody

Cat No: HR1AP4854

For research use only

Overview

Product Name	ZNF771 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GenelD?Human?	ZNF771
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF771; Zinc finger protein 771; Mesenchymal stem cell protein DSC43
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 8 C2H2-type zinc fingers.,