



# ZNF596 Polyclonal Antibody

Cat No: HR1AP6105

For research use only

## Overview

|                       |                                                                                                                                                                                                                     |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | ZNF596 Polyclonal Antibody                                                                                                                                                                                          |
| Source                | Rabbit                                                                                                                                                                                                              |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                      |
| Species Reactivity    | Human                                                                                                                                                                                                               |
| Recommended Dilutions |                                                                                                                                                                                                                     |
| Immunogen             |                                                                                                                                                                                                                     |
| Species               | Rabbit                                                                                                                                                                                                              |
| Storage               | -20°C/1 year                                                                                                                                                                                                        |
| Isotype               |                                                                                                                                                                                                                     |
| Clonality             |                                                                                                                                                                                                                     |
| Concentration         | 1 mg/ml                                                                                                                                                                                                             |
| Observed band         | 58kDa                                                                                                                                                                                                               |
| GeneID?Human?         | ZNF596                                                                                                                                                                                                              |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                     |
| Cellular localization |                                                                                                                                                                                                                     |
| Alternative Names     | ZNF596; Zinc finger protein 596                                                                                                                                                                                     |
| Background            | function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 11 C2H2-type zinc fingers., |