

CD68 Mouse mAb

Cat No: HR1AM2273

For research use only

Overview

Product Name	CD68 Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	968
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD68 molecule, GP110, Macrosialin, SCARD1
Background	The CD68 antigen is a 37kD transmembrane protein that is post-translationally glycosylated to give a protein of 87-115kD. CD68 is specifically expressed by tissue macrophages, Langerhans cells and at low levels by dendritic cells. It could play a role in phagocytic activities of tissue macrophages, both in intracellular lysosomal metabolism and extracellular cell-cell and cell-pathogen interactions. It binds to tissue- and organ-specific lectins or selectins, allowing homing of macrophage subsets to particular sites. Rapid recirculation of CD68 from endosomes and lysosomes to the plasma membrane may allow macrophages to crawl over selectin bearing substrates or other cells.