



ZNF337 Polyclonal Antibody

Cat No: HR1AP4847

For research use only

Overview

Product Name	ZNF337 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	86kDa
GeneID?Human?	ZNF337
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF337; Zinc finger protein 337
Background	zinc finger protein 337(ZNF337) Homo sapiens This gene encodes a zinc finger domain containing protein. The function of this protein has yet to be determined. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014],