



ZNF225 Polyclonal Antibody

Cat No: HR1AP4845

For research use only

Overview

Product Name	ZNF225 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	82kDa
GeneID?Human?	ZNF225
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF225; Zinc finger protein 225
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 18 C2H2-type zinc fingers.,