



## ZNF134 Polyclonal Antibody

Cat No: HR1AP4844

For research use only

### Overview

Product Name	ZNF134 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GenelD?Human?	ZNF134
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF134; Zinc finger protein 134
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 11 C2H2-type zinc fingers.,