

Lactoferrin Mouse mAb

Cat No: HR1AM2272

For research use only

Overview

Product Name	Lactoferrin Mouse mAb
Source	Mouse
Applications	ELISA
Species Reactivity	Human
Recommended Dilutions	ELISA: 1:5,000-10,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	77kDakDa
GeneID?Human?	4057
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CKRX antibody, GIG12 antibody, HLF2 antibody, LF antibody, LFT antibody, Talalactoferrin antibody
Background	Lactoferrin (LF), also known as lactotransferrin (LTF), is a multifunctional protein of the transferrin family. Lactoferrin is one of the components of the immune system of the body; it has antimicrobial activity (bacteriocide, fungicide) and is part of the innate defense.