

ZBTB40 Polyclonal Antibody

Cat No: HR1AP4836

For research use only

Overview

Product Name	ZBTB40 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	140kDa
GeneID?Human?	ZBTB40
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZBTB40; KIAA0478; Zinc finger and BTB domain-containing protein 40
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 12 C2H2-type zinc fingers.,