



ZC3H13 Polyclonal Antibody

Cat No: HR1AP4838

For research use only

Overview

Product Name	ZC3H13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	200kDa
GenID?Human?	ZC3H13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZC3H13; KIAA0853; Zinc finger CCCH domain-containing protein 13
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 C3H1-type zinc finger.,



bioelsa