

XRN2 Polyclonal Antibody

Cat No: HR1AP4828

For research use only

Overview

Product Name	XRN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	108kDa
GeneID?Human?	XRN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	XRN2; 5'-3' exoribonuclease 2; DHM1-like protein; DHP protein
Background	5'-3' exoribonuclease 2(XRN2) Homo sapiens This gene encodes a 5''-3'' exonuclease that promotes transcription termination at cotranscriptional cleavage sites. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015],