



TBX10 Polyclonal Antibody

Cat No: HR1AP4677

For research use only

Overview

Product Name	TBX10 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GenelD?Human?	TBX10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TBX10; TBX7; T-box transcription factor TBX10; T-box protein 10
Background	T-box 10(TBX10) <i>Homo sapiens</i> This gene encodes a member of the T-box family of transcription factors. These transcription factors share a DNA-binding domain called the T-box, and play a role in several developmental processes including early embryonic cell fate and organogenesis. The encoded protein is a member of the T-box 1 subfamily. Mutations in this gene are thought to be a cause of isolated cleft lip with or without cleft palate. [provided by RefSeq, Nov 2010],