




Spi-C Polyclonal Antibody

Cat No: HR1AP4595

For research use only

Overview

Product Name	Spi-C Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GenID?Human?	SPIC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SPIC; Transcription factor Spi-C
Background	Spi-C transcription factor(SPIC) Homo sapiens The protein encoded by this gene regulates the development of red pulp macrophages, which are necessary for iron homeostasis and the recycling of red blood cells. [provided by RefSeq, Aug 2016],



bioelsa