



# Spi-C Polyclonal Antibody

Cat No: HR1AP4595

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                  |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Spi-C Polyclonal Antibody                                                                                                                                                                                                                        |
| Source                | Rabbit                                                                                                                                                                                                                                           |
| Applications          | WB,IF,ELISA                                                                                                                                                                                                                                      |
| Species Reactivity    | Human                                                                                                                                                                                                                                            |
| Recommended Dilutions |                                                                                                                                                                                                                                                  |
| Immunogen             |                                                                                                                                                                                                                                                  |
| Species               | Rabbit                                                                                                                                                                                                                                           |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                     |
| Isotype               |                                                                                                                                                                                                                                                  |
| Clonality             |                                                                                                                                                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                          |
| Observed band         | 29kDa                                                                                                                                                                                                                                            |
| GeneID?Human?         | SPIC                                                                                                                                                                                                                                             |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                  |
| Cellular localization |                                                                                                                                                                                                                                                  |
| Alternative Names     | SPIC; Transcription factor Spi-C                                                                                                                                                                                                                 |
| Background            | Spi-C transcription factor(SPIC) Homo sapiens The protein encoded by this gene regulates the development of red pulp macrophages, which are necessary for iron homeostasis and the recycling of red blood cells. [provided by RefSeq, Aug 2016], |