

SCAND1 Polyclonal Antibody

Cat No: HR1AP4521

For research use only

Overview

Product Name	SCAND1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	SCAND1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SCAND1; SDP1; SCAN domain-containing protein 1
Background	SCAN domain containing 1(SCAND1) Homo sapiens This gene encodes a SCAN box domain-containing protein. The SCAN domain is a highly conserved, leucine-rich motif of approximately 60 aa originally found within a subfamily of zinc finger proteins. This gene belongs to a family of genes that encode an isolated SCAN domain, but no zinc finger motif. This protein binds to and may regulate the function of the transcription factor myeloid zinc finger 1B. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Jan 2011],