

Flotillin-1 Mouse mAb

Cat No: HR1AM2268

For research use only

Overview

Product Name	Flotillin-1 Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	49kDakDa
GeneID?Human?	10211
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FLOT 1 antibody, Integral membrane component of caveolae antibody, Reggie 2 antibody
Background	Two flotillin members have been described, flotillin-1 and flotillin-2. May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.