

RoXaN Polyclonal Antibody

Cat No: HR1AP4502

For research use only

Overview

Product Name	RoXaN Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	111kDa
GeneID?Human?	ZC3H7B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZC3H7B; KIAA1031; Zinc finger CCCH domain-containing protein 7B; Rotavirus 'X'-associated non-structural protein; RoXaN
Background	zinc finger CCCH-type containing 7B(ZC3H7B) Homo sapiens This gene encodes a protein that contains a tetratricopeptide repeat domain. The encoded protein also interacts with the rotavirus non-structural protein NSP3. [provided by RefSeq, Jul 2008],