



RFX3 Polyclonal Antibody

Cat No: HR1AP4456

For research use only

Overview

Product Name	RFX3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GenID?Human?	RFX3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RFX3; Transcription factor RFX3; Regulatory factor X 3
Background	regulatory factor X3(RFX3) <i>Homo sapiens</i> This gene is a member of the regulatory factor X gene family, which encodes transcription factors that contain a highly-conserved winged helix DNA binding domain. The protein encoded by this gene is structurally related to regulatory factors X1, X2, X4, and X5. It is a transcriptional activator that can bind DNA as a monomer or as a heterodimer with other RFX family members. Multiple transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2013].