

PRDM12 Polyclonal Antibody

Cat No: HR1AP4369

For research use only

Overview

Product Name	PRDM12 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	PRDM12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PRDM12; PFM9; PR domain zinc finger protein 12; PR domain-containing protein 12
Background	function:May be involved in transcriptional regulation.,similarity:Contains 1 SET domain.,similarity:Contains 3 C2H2-type zinc fingers.,