



OSR2 Polyclonal Antibody

Cat No: HR1AP5770

For research use only

Overview

Product Name	OSR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	OSR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OSR2; Protein odd-skipped-related 2
Background	odd-skipped related transcription factor 2(OSR2) Homo sapiens OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).[supplied by OMIM, Mar 2008],