



## NRIP3 Polyclonal Antibody

Cat No: HR1AP4098

For research use only

### Overview

Product Name	NRIP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GenID?Human?	NRIP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NRIP3; C11orf14; Nuclear receptor-interacting protein 3; Sarcoma antigen NY-SAR-105
Background	NRIP3 (Nuclear Receptor Interacting Protein 3) is a Protein Coding gene. Gene Ontology (GO) annotations related to NRIP3 include aspartic-type endopeptidase activity. An important paralog of NRIP3 is NRIP2.