



Nrl Polyclonal Antibody

Cat No: HR1AP4099

For research use only

Overview

Product Name	Nrl Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GenID?Human?	NRL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NRL; D14S46E; Neural retina-specific leucine zipper protein; NRL
Background	neural retina leucine zipper(NRL) <i>Homo sapiens</i> This gene encodes a basic motif-leucine zipper transcription factor of the Maf subfamily. The encoded protein is conserved among vertebrates and is a critical intrinsic regulator of photoreceptor development and function. Mutations in this gene have been associated with retinitis pigmentosa and retinal degenerative diseases. [provided by RefSeq, Jul 2008],