



## MN1 Polyclonal Antibody

Cat No: HR1AP3918

For research use only

### Overview

Product Name	MN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	135kDa
GenID?Human?	MN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MN1; Probable tumor suppressor protein MN1
Background	MN1 proto-oncogene, transcriptional regulator(MN1) Homo sapiens Meningioma 1 (MN1) contains two sets of CAG repeats. It is disrupted by a balanced translocation (4;22) in a meningioma, and its inactivation may contribute to meningioma 32 pathogenesis. [provided by RefSeq, Jul 2008],