



SRT-Tag Mouse mAb

Cat No: HR1AM2265

For research use only

Overview

Product Name	SRT-Tag Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	N/A
Recommended Dilutions	WB 1:5,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SRT Tag antibody
Background	The SRT-tag is a peptide epitope comprised of 10 amino acids(TFIGAIATDT) derived from the crystalline surface layer protein of Rickettsia typhi1.