

Med18 Polyclonal Antibody

Cat No: HR1AP3876

For research use only

Overview

Product Name	Med18 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	MED18
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MED18; Mediator of RNA polymerase II transcription subunit 18; Mediator complex subunit 18; p28b
Background	mediator complex subunit 18(MED18) Homo sapiens MED18 is a component of the Mediator complex, which is a coactivator for DNA-binding factors that activate transcription via RNA polymerase II (Sato et al., 2003 [PubMed 12584197]).[supplied by OMIM, Oct 2008],