



Med18 Polyclonal Antibody

Cat No: HR1AP3876

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | Med18 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 24kDa |
| GenID?Human? | MED18 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | MED18; Mediator of RNA polymerase II transcription subunit 18; Mediator complex subunit 18; p28b |
| Background | mediator complex subunit 18(MED18) Homo sapiens MED18 is a component of the Mediator complex, which is a coactivator for DNA-binding factors that activate transcription via RNA polymerase II (Sato et al., 2003 [PubMed 12584197]) [supplied by OMIM, Oct 2008], |