

? II tubulin Mouse mAb Cat No: HR1AM2262

For research use only

Overview

Product Name	? II tubulin Mouse mAb
Source	Mouse
Applications	WB, IP
Species Reactivity	Rat,Mouse,Human
Recommended Dilutions	WB 1:50,000-100,000 IP 1:200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgGq
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	50kDakDa
GeneID?Human?	7280
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	class II beta tubulin isotype antibody, TUBB 2 antibody, Tubulin beta 2 antibody
Background	Tubulin is the major constituent of microtubules. It exists as a heterodimer consiting of an alpha and a beta subunit. Class II beta tubulin is the major form of tubulin beta in neurons, although it is not neuronal specific. It is found in many other tissues including lung tissue and Schwann cells.