




Islet-2 Polyclonal Antibody

Cat No: HR1AP3762

For research use only

Overview

Product Name	Islet-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GenID?Human?	ISL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ISL2; Insulin gene enhancer protein ISL-2; Islet-2
Background	function:Transcriptional factor that defines subclasses of motoneurons that segregate into columns in the spinal cord and select distinct axon pathways.,similarity:Contains 1 homeobox DNA-binding domain.,similarity:Contains 2 LIM zinc-binding domains.,



bioelsa