



Id4 Polyclonal Antibody

Cat No: HR1AP3702

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | Id4 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 24kDa |
| GenelD?Human? | ID4 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | ID4; BHLHB27; DNA-binding protein inhibitor ID-4; Class B basic helix-loop-helix protein 27; bHLHb27; Inhibitor of DNA binding 4 |
| Background | inhibitor of DNA binding 4, HLH protein(ID4) Homo sapiens This gene encodes a member of the inhibitor of DNA binding (ID) protein family. These proteins are basic helix-loop-helix transcription factors which can act as tumor suppressors but lack DNA binding activity. Consequently, the activity of the encoded protein depends on the protein binding partner. [provided by RefSeq, Dec 2011], |