

HSFX1 Polyclonal Antibody

Cat No: HR1AP3684

For research use only

Overview

Product Name	HSFX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	HSFX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HSFX1; LW-1; HSFX2; Heat shock transcription factor; X-linked
Background	similarity:Belongs to the HSF family.,tissue specificity:Testis-specific.,