



HKR1 Polyclonal Antibody

Cat No: HR1AP3646

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | HKR1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,IF,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 73kDa |
| GenelD?Human? | HKR1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | HKR1; ZNF875; Krueppel-related zinc finger protein 1; Protein HKR1; Zinc finger protein 875 |
| Background | function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 13 C2H2-type zinc fingers., |