



FoxN2 Polyclonal Antibody

Cat No: HR1AP6908

For research use only

Overview

Product Name	FoxN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	FOXN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXN2; HTLF; Forkhead box protein N2; Human T-cell leukemia virus enhancer factor
Background	forkhead box N2(FOXN2) Homo sapiens This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat. [provided by RefSeq, Jul 2008],