

FoxB2 Polyclonal Antibody

Cat No: HR1AP3472

For research use only

Overview

Product Name	FoxB2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	FOXB2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXB2; Forkhead box protein B2
Background	function:Transcription factor.,similarity:Contains 1 fork-head DNA-binding domain.,