

FoxD4L1 Polyclonal Antibody

Cat No: HR1AP6211

For research use only

Overview

Product Name	FoxD4L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	FOXD4L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXD4L1; Forkhead box protein D4-like 1; FOXD4-like 1
Background	forkhead box D4-like 1(FOXD4L1) Homo sapiens This gene is a member of the forkhead/winged-helix (FOX) family of transcription factors with highly conserved FOX DNA-binding domains. Members of the FOX family of transcription factors are regulators of embryogenesis and may play a role in human cancer. This gene lies in a region of chromosome 2 that surrounds the site where two ancestral chromosomes fused to form human chromosome 2. This region is duplicated elsewhere in the human genome, primarily in subtelomeric and pericentromeric locations, thus mutiple copies of this gene have been found. [provided by RefSeq, Jul 2008],