




## Fliz1 Polyclonal Antibody

Cat No: HR1AP3461

For research use only

### Overview

Product Name	Fliz1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GenID?Human?	ZC3H8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZC3H8; ZC3HDC8; Zinc finger CCCH domain-containing protein 8
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 3 C3H1-type zinc fingers.,



**bioelsa**