




FoxB1/2 Polyclonal Antibody

Cat No: HR1AP3471

For research use only

Overview

Product Name	FoxB1/2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GenID?Human?	FOXB1/FOXB2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXB1; FKH5; Forkhead box protein B1; Transcription factor FKH-5; FOXB2; Forkhead box protein B2
Background	similarity:Contains 1 fork-head DNA-binding domain.,



bioelsa