



FBP3 Polyclonal Antibody

Cat No: HR1AP3446

For research use only

Overview

Product Name	FBP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GenID?Human?	FUBP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FUBP3; FBP3; Far upstream element-binding protein 3; FUSE-binding protein 3
Background	function:May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 4 KH domains.,tissue specificity:Detected in a number of cell lines.,