



# FBP3 Polyclonal Antibody

Cat No: HR1AP3446

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                          |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | FBP3 Polyclonal Antibody                                                                                                                                                                                                                                                 |
| Source                | Rabbit                                                                                                                                                                                                                                                                   |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                           |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                          |
| Recommended Dilutions |                                                                                                                                                                                                                                                                          |
| Immunogen             |                                                                                                                                                                                                                                                                          |
| Species               | Rabbit                                                                                                                                                                                                                                                                   |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                             |
| Isotype               |                                                                                                                                                                                                                                                                          |
| Clonality             |                                                                                                                                                                                                                                                                          |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                  |
| Observed band         | 60kDa                                                                                                                                                                                                                                                                    |
| GeneID?Human?         | FUBP3                                                                                                                                                                                                                                                                    |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                          |
| Cellular localization |                                                                                                                                                                                                                                                                          |
| Alternative Names     | FUBP3; FBP3; Far upstream element-binding protein 3; FUSE-binding protein 3                                                                                                                                                                                              |
| Background            | function:May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 4 KH domains.,tissue specificity:Detected in a number of cell lines., |