



# Elf-1 Polyclonal Antibody

Cat No: HR1AP6199

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                           |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Elf-1 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                 |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                    |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                                                                                                                                                                                                                         |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                                                                                     |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                           |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                           |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                    |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                              |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                           |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                           |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                   |
| Observed band         | 67kDa                                                                                                                                                                                                                                                                                                                                                                                     |
| GeneID?Human?         | ELF1                                                                                                                                                                                                                                                                                                                                                                                      |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                           |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                           |
| Alternative Names     | ELF1; ETS-related transcription factor Elf-1; E74-like factor 1                                                                                                                                                                                                                                                                                                                           |
| Background            | E74 like ETS transcription factor 1(ELF1) Homo sapiens This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2009], |