



Elf-1 Polyclonal Antibody

Cat No: HR1AP6199

For research use only

Overview

Product Name	Elf-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GenID?Human?	ELF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ELF1; ETS-related transcription factor Elf-1; E74-like factor 1
Background	E74 like ETS transcription factor 1(ELF1) Homo sapiens This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2009],