



## ECP Polyclonal Antibody

Cat No: HR1AP3332

For research use only

### Overview

Product Name	ECP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GenID?Human?	RNASE3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RNASE3; ECP; RNS3; Eosinophil cationic protein; ECP; Ribonuclease 3; RNase 3
Background	ribonuclease A family member 3(RNASE3) Homo sapiens The protein encoded by this gene belongs to the pancreatic ribonuclease family, a subset of the ribonuclease A superfamily. The protein exhibits antimicrobial activity against pathogenic bacteria [provided by RefSeq, Oct 2014],