

DQX1 Polyclonal Antibody

Cat No: HR1AP3308

For research use only

Overview

Product Name	DQX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	DQX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DQX1; ATP-dependent RNA helicase DQX1; DEAQ box polypeptide 1
Background	similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,