



DNA pol ? Polyclonal Antibody

Cat No: HR1AP3295

For research use only

Overview

Product Name	DNA pol ? Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	165kDa
GenelD?Human?	POLA1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	POLA1; POLA; DNA polymerase alpha catalytic subunit; DNA polymerase alpha catalytic subunit p180
Background	DNA polymerase alpha 1, catalytic subunit(POLA1) <i>Homo sapiens</i> This gene encodes the catalytic subunit of DNA polymerase, which together with a regulatory and two primase subunits, forms the DNA polymerase alpha complex. The catalytic subunit plays an essential role in the initiation of DNA replication. [provided by RefSeq, Mar 2010],