

# DDX52 Polyclonal Antibody

Cat No: HR1AP3268

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | DDX52 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                      |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Observed band         | 67kDa                                                                                                                                                                                                                                                                                                                                                                                                                          |
| GeneID?Human?         | DDX52                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Alternative Names     | DDX52; ROK1; HUSSY-19; Probable ATP-dependent RNA helicase DDX52; ATP-dependent RNA helicase ROK1-like; DEAD box protein 52                                                                                                                                                                                                                                                                                                    |
| Background            | caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Several sequencing problems.,similarity:Belongs to the DEAD box helicase family. DDX52/ROK1 subfamily.,similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain., |