



DDX52 Polyclonal Antibody

Cat No: HR1AP3268

For research use only

Overview

Product Name	DDX52 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GenID?Human?	DDX52
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DDX52; ROK1; HUSSY-19; Probable ATP-dependent RNA helicase DDX52; ATP-dependent RNA helicase ROK1-like; DEAD box protein 52
Background	caution: The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., sequence caution: Several sequencing problems., similarity: Belongs to the DEAD box helicase family. DDX52/ROK1 subfamily., similarity: Contains 1 helicase ATP-binding domain., similarity: Contains 1 helicase C-terminal domain.,