

DDX51 Polyclonal Antibody

Cat No: HR1AP3267

For research use only

Overview

Product Name	DDX51 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	72kDa
GeneID?Human?	DDX51
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DDX51; ATP-dependent RNA helicase DDX51; DEAD box protein 51
Background	domain:The Q motif is unique to and characteristic of the DEAD box family of RNA helicases and controls ATP binding and hydrolysis.,function:ATP-binding RNA helicase involved in the biogenesis of 60S ribosomal subunits.,similarity:Belongs to the DEAD box helicase family. DDX51/DBP6 subfamily.,similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,